

**World Firstly Developed
Extracorporeal non invasive active Biofeed back
device to do Kegel Exercise for urinary
incontinence, fecal incontinence patients
HnJ-7000 for Hospitals, Clinics**



FURUN HEALTHCARE

Composition of the product



Main body



Abdominal pressure
detection Belt



Data cable



Operating computer
(optional)



Printer(optional)

- Supplying package per unit

1. Main Body (Chair) ... [1EA]

2. Abdominal Pressure detection Belt [1EA]

3. Data cable [1EA]

4. computer (optional) [1EA]

5. Printer (optional) [1EA]




Product spec

● Pelvic Floor Muscles Training System (Model : HnJ-7000)

Division	Contents	Remarks
Input power	AC 100 ~ 240V, 50Hz / 60Hz	
Power Consumption	Max 200W	
Range of strength detection	0 ~ 25kgf	
Size	640× 820× 1,100mm (W-D-H)	
Weight	53kg	
Display	Control PC, 4GHz	Optional

Comparison table with global competitors

- For Hospital Use

Classification	Biofeedback device (Probe insertion type)	Electromagnetic device	Biofeedback device (without probe insertion)
	Global Competition Company		Unique product in the world
Shape & Function			
How it works	Insert probe into vagina / Active Exercise	No insertion / Magnetic stimulation	No insertion / Active Exercise
Ease of use	Take off / Probe insertion (X)	Put on / No insertion (O)	Put on / No insertion (O)
Muscle strengthening	No-load training (O) (active type)	Magnetic stimulation (X) (passive type)	Excellent (O) weight training (active)
Muscles control ability	good (Δ)	None (X)	Excellent (O)
Learning	Difficulty (Δ)	None (X)	Easy to learn (O)

Features

Convenience

Realized exactly same **pelvic floor muscle exercise with cloths on** which was possible by inserting a probe into the female vagina in the past

Effective

Weight training for your pelvic muscles similar to arm workouts with dumbbell

Motivation

Real-time monitoring of pelvic floor muscles contraction and relaxation

Various treatment

Urinary Incontinence , Feces , Constipation, Irritable bladder
Pelvic Floor Muscle disorder

Anyone, anywhere

Convenient for young and old, and male and female