

World Firstly Developed
Extracorporeal non invasive active Biofeed back
device to do Kegel Exercise
HnJ-350E for Home ,Personal use



FURUN HEALTHCARE

Composition of the product



Main body



Smart phone or Tablet PC
(Need to use customer's smartphone)



Charging adaptor



Data cable



User manual /Quick manual



Wet tissue

● Supplying package per unit



- | | |
|---|---------------------------------|
| 1. Main Body (Cushion) ... [1EA] | 2. Charging Adaptor [1EA] |
| 3. Data Cable(Main Body <---> Smart phone, 1.5m) | [1EA] |
| 4. User/Quick Manual [1EA/Each] | 5. Wet Tissue [1EA] |

Product spec

Pelvic Floor Muscles Training System (Model : HnJ-350)




| | | |
|------------------------------------|-------------------------------|--------------------------|
| 1. Main Body (Cushion) | Input Power | DC 5V |
| | Battery Capacity | 1,800mAh |
| | measurement range of strength | 0 ~ 20kgf |
| | Dimensions | 400 x 300 x 50mm (W-D-H) |
| | Weight | 1.3kg |
| 2. Display part | Smartphone | Android Phone / i-phone |
| | Tablet PC | Option : |
| 3 Operation environment | Temperature | 10°C - 50°C |
| | Humidity | 0% - 75% RH |
| 4. Transport and storage condition | Temperature | -20°C - 70°C |
| | Humidity | 0% - 75% RH |
| | Air Pressure | 700 – 1060 mbar |

Comparison table with global competitors


| | Furun Healthcare | | Global competitors | |
|-----------------|---|---|--|---|
| | Bio feedback device with cloths on (World firstly developed) | | Bio feed back device, Insertion probe into the vagina (most widely using now) | |
| Product type |  | |  <p>Kegel Exercises: Contract your pelvic floor muscles for three seconds, then relax the muscles for three seconds. Do this 10-15 times several times a day. Although shown here while lying down, these exercises can be done during a variety of daily activities, such as sitting in a meeting, while stopped in your car at a traffic light or when talking on the phone.</p> | |
| Exercise scheme | Pelvic floor muscle training and low frequency stimulation together with cloths on | | Pelvic floor muscle training with inserting probe into the vagina (Graph type) | |
| Convenience | Can use with cloths on / Non insertion type | 0 | Take off the cloths / probe insertion type | X |
| Motivation | Various kinds of entertainments(Games & Fun) | 0 | Simple graph/ number | X |
| Treatment way | Muscle training (Kegel exercise) Low frequency stimulation (urinary incontinence+ reduce pain) | 0 | Simple muscle training | △ |
| Target | For Man, Woman, Old, Youth every body, Anywhere | 0 | Woman, only separated space | X |

Function of the product

- Kegels exercise mode

| | | | |
|--------------------|---|---|--|
| <p>screen</p> |  |  |  |
| | <p>Graph mode</p> | <p>Game 1</p> | <p>Game 2</p> |
| <p>Description</p> | <p>Graph mode appearing various short and long twitch exercise graph.</p> | <p>Designed to score coins by exercise pelvic floor muscles to improve pelvic floor muscle strength</p> | <p>Improving controlling ability for the pelvic floor muscles with contraction and relaxation.</p> |

- Low-frequency stimulation mode



Low frequency screen



Pelvic stimulation



Pain Relief

Features

Same principle
with hospital
device

The Same bio feed back device principle as hospital device which are using in big hospitals, clinics now

Convenience

Realized exactly same pelvic floor muscle exercise with cloths on which was possible by inserting a probe into the female vagina in the past

Various
functions

Exercise therapy through Kegel exercise,
Nerve stimulation therapy through low frequency stimulation.
waist, shoulder, knee ,body, etc.. low frequency pain treatment

Easy to
use

Mother, father, grand mother, grand father, whoever...
Everybody can easily use.

Maximum
Treatment
effectiveness

Maximizing treatment effect by customized
exercise therapy and stimulation therapy to user.